

Legacy Cyclone Taper

Nickel Titanium Rotary Files

RX Only - For Dental Use Only




Single Use Only



Consult Instructions For Use



Sterilize Before Use



Autoclave Temperature and Time

Composition

Operative portion is composed of nickel titanium alloy.

Indications for use

For the removal of dentin and shaping of the root canal.

Warnings

Files are provided non-sterile, and must be sterilized before use. This product contains nickel and should not be used for individuals with known allergic sensitivity to this material.

Contraindications

None Known

Adverse Reactions

None Known

Precautions

1. As with all new products, you must exercise caution until you become proficient in its use.
2. A slow-speed hand piece is required for rotary file use.
3. Operate the hand piece at 500 RPM (Revolutions per minute) with a torque setting of 300 g-cm.
4. Straight-line access is a prerequisite for proper endodontic treatment.
5. Always utilize minimal apical pressure. Do not force the files down the canal.
6. Clean file flutes after each insertion into the root canal.
7. Frequently irrigate and lubricate the canal throughout the procedure.
8. Take each rotary file to length only one time and for no more than one second.
9. Exercise caution in the apical area and around significant curvatures.
10. Rotary files are single patient use devices.
11. If the file is too large and taken to length it will increase the risk of canal transportation and file separation.

Sterilization

Files must be sterilized before use. ANSI/ADA Specification 28 recommends:

- Scrub the instruments with soap and warm water
- Rinse thoroughly with distilled or deionized water
- Allow to air dry
- Place the instruments, unwrapped, in the autoclave tray
- Use fresh distilled or deionized water
- Steam Autoclave at 136°C (+/- 2°C) for 20 minutes
- Rotary files are single patient use instruments
- Recommended file disposal: Place used files in Biohazard Sharps container

Step-by-Step IFU

Straight-Line Access

1. Create straight line access to canal orifices.
2. Locate canals and explore using small stainless steel hand instrument. Minimum size #15 K-File required prior to rotary use. Apex Access Kit rotary files size 15 and 20 can also be used after #10 K-File in lieu of additional hand instruments.

Crown down Shaping and Cleaning 04 Taper

For smaller canals start with a 30/04 file. Take the 30/04 to resistance or working length. If resistance is encountered before working length is obtained, go to the next smaller instrument following the same protocol until working length is achieved. Between each rotary file recapitulate with a #10 or #15 tip hand file to maintain glide path and help irrigate to the canal terminus. For larger canals start with a 40/04 and proceed using the Crown Down technique to resistance or working length is achieved.

The same technique applies for all 06 tapered instruments if required.