



# Legacy 7



The Legacy 7 is a series of 7 files (17, 20, 25, 30, 35, 40, and 45 tip sizes) available in .04 and .06 tapers. The patented Controlled Memory process developed by D&S Dental allows the Legacy 7 series to reach maximum flexibility and ensure a safe and effective use. This series is available in 21mm, 25mm, and 31mm working lengths. These instruments are used with familiar techniques and sequences that are comparable to other products on the market.

**Compare to:** EdgeFile-X7™ (Made in France)



Cross Section



855-QDP-0220 | [www.QDPendo.com](http://www.QDPendo.com)

