


Legacy Cyclone Taper
Legacy Cyclone Taper CM



The Legacy Cyclone Taper instruments are designed with a triangular cross section and variable pitch that provides optimal cutting efficiency and flexibility. The Cyclone Taper is offered in both .04 and .06 tapers. Available working lengths are 21mm, 25mm, and 31mm. Cyclone Taper files utilize a Nickel plated rotary handle which are manufactured to ISO specifications. All files have non cutting tips to prevent ledging. This system is available in Standard NiTi and Controlled Memory to increase flexibility.

Compare to: Vortex® and Vortex Blue® (Made in Germany and Switzerland)



Cross Section

