


Legacy MT & Legacy MT CM



The Legacy MT is a Multi Tapered instrument designed to provide a flexible, safe, and efficient cutting combination. The MT file system includes 3 shaping files (SX, S1, S2) and 5 Finishing files (F1, F2, F3, F4, F5). The MT is available in 21mm, 25mm, and 31mm working lengths with the exception of the SX file which has a 19mm working length. All files have non cutting tips to prevent ledging. This system is available in Standard NiTi and Controlled Memory to increase flexibility.

Compare to: ProTaper® and ProTaper Gold® (*Made in Germany and Switzerland*)
EdgeTaper™ and EdgeTaper Platinum™ (*Made in France*)



Cross Section



855-QDP-0220 | www.QDPendo.com

